



Jewell Nursing Solutions

Pressure Injury Care & Prevention

The FootBed Heel Floatation Cushion

Pressure Injury Relief

**Taking the pressure off patients
and caregivers since 2012**

Finally - Real, practical, and effective solutions for pressure injury prevention and treatment everyone can use.



**CLICK
HERE**
to view
video



**Jewell
Nursing
Solutions**

(650) 294-8557 | support@jewellnursingsolutions.com | www.jewellnursingsolutions.com
Bedsore Rescue™ is patented and trademarked by Jewell Nursing Solutions. All rights reserved.

Footbed Heel Floater

Specifications

Overall Dimensions	25" x 19" x 5.5"		
Weight	3lb.		
Inner Cushion Material	<ul style="list-style-type: none"> • Center Foam Core: High resolution, medium density, polyurethane foam • Outer Rail: Chopped Microfiber for longer lasting fill volume <p>Latex free, Hypoallergenic, Contains no metals, no zippers</p>		
Cover	<ul style="list-style-type: none"> • Top and Side Rails: 100% ChemflexII® Polyester knit 'gurney pad' fabric • Bottom: Roughened Solid Polyurethane • Moisture tolerance: Proof Resistant to 150lbs/per sq. inch • Moisture vapor exchange: Breathability-194 G/m2/24hr (soft feel, skin safe) • Antimicrobial Treatment: Ultra Fresche formula, hypoallergenic • Construction: Fine stitch fixed. Contains no metals, no zippers, and no seams where skin touches cushion 		
Disinfection	<p>This product is reusable. Anti-microbially treated to inhibit colonization. Wipe or rinse to clean and disinfect between uses. Compatible with bleach, hydrogen peroxide, chlorhexidine and alcohol based disinfectant wipes, or soap & water. Do not machine wash. Do not submerge. Discard if cover is damaged or contamination is suspected. We recommend discarding after use in COVID-19 isolation. Do not clean with caustic cleaners made for degradation or stripping.</p>		
	Product Name Footbed Heel Floater	SKU FB200-HF	UPC 661799582643
Ordering	<p>For orders less than six cushions, go to: www.jewellnursingsolutions.com For contract and bulk orders, send PO to: support@jewellnursingsolutions.com Support phone number is (650) 294-8557</p>		
Guarantee & Warranty	<p>Three months warranty for materials & construction under normal wear and tear (Excludes damage from sharp objects, corrosive fluids, submersion in fluids (laundering, tubs), and exposure to temperatures >200 degrees). 10-day satisfaction guarantee or your money back when purchasing two or less cushions.</p>		

Meets flammability requirements of electronic & appliance repair, home furnishing, and thermal insulation technical bulletin 117-2013. Contains NO added flame-retardant chemicals.

Jewell Nursing Solutions is a nurse owned and operated company dedicated to bringing effective pressure injury prevention products and services to the entire caregiving community. Bedsore Rescue® Cushion (formerly The BackBone™) is proudly designed, patented, and trademarked by Jewell Nursing Solutions. A certified woman owned business celebrating 10 years of helping people stop pressure injuries. Made in the USA by certified manufacturer reg #CA39930 (IL). Distributed by Jewell Nursing Solutions

How to Use the Footbed Heel Floatation Cushion

1. Side rail for support. Gap space between rail and center pad for floating heels and ankle bones (aka malleolous bony prominences)



2. Low profile design provides just enough lift to relieve pressure but not dangle ankles for ultimate comfort. Does not affect venous and arterial blood flow. May place pillows underneath cushion when leg elevation is desired. Ankle support will still be achieved.



3. Allow foot to rest so that both heel and ankle bones are floating, but that there is gentle support under achillies and side of foot



4. Float heels over gaps while supporting bottom of feet.

Okay to put small pillow between feet for more support.



5. When legs spread to outside of cushion, ankles and side of leg are supported in corners while heels float over gap spaces.



6. Rotate cushion to accommodate natural leg positions. Float heels & ankle bones and support feet without elevating legs.



7. Rotate cushion so that rails support underneath knees, and heels and ankle bony prominences are floating. It's okay to place pillow between knees for comfort.



COMFORT TIPS

Everybody is different. To find a position that can achieve maximum comfort & best possible pressure relief, take just a few seconds each time you are in the room to adjust/ rotate the cushion to accomodate floatation and comfort.

NURSING ADVICE

The footbed is designed to float the heels and support the ankles without affecting blood flow. Generally a low profile support is preferred for wound healing. But follow your doctor or woundcare nurses recommendations if the patient has swelling or venous stasis syndrome. We recommend alternating between an elevated position and a level position. If you are having difficulty with turning and repositioning, visit our website for education and training on how to turn and reposition for best possible pressure relief and comfort. Or feel free to contact us directly if you need specific nursing advice or want to learn about a few tricks of the trade!

www.jewellnursingsolutions.com
support@jewellnursingsolutions.com
650-297-8557