



The BedSore Rescue™ Positioning Wedge

Selection Guide and Specifications



	<p>Standard Size 27" L x 12" W Wedge 5" at back to 2" at front tips</p>	<p>Large Size 28" L x 12" W Wedge 6" at back to 2.5" at front tips</p>	<p>Bariatric Size 29" L x 18" W Wedge 8" at back to 3" at front tips</p>																					
Persons Weight	80 to 220lbs	221-350 pounds	351- 600 pounds																					
Persons Height	5'0" to 6'.5"																							
Material	<p>Inner Core: High resolution medium density polyurethane foam Support layer: Pressure distributive high density memory foam Latex Free, Contains no metals, no zippers</p>																							
Fabric Cover	<p>Home Version: 100% ChemflexII Polyester knit "gurney pad" fabric Non-skid general use: Top: Chemflex II ® Polyurethane coated knit "gurney pad" fabric Non-skid Advanced: Top: Carrflex CFX71EF® Polyurethane coated interlock polyester knit Bottom (Both general use and advance): Roughened Solid Polyurethane</p> <table border="1" data-bbox="396 945 1507 1134"> <thead> <tr> <th></th> <th>Moisture Resistance</th> <th>Vapor Exchange</th> <th>Stretch Way</th> <th>Tear/Break Strength</th> <th>Seams Type</th> <th>Antimicrobial Treatment</th> </tr> </thead> <tbody> <tr> <td>ChemflexII:</td> <td>150lbs/sq.in</td> <td>194 G/m2/24hr</td> <td>2</td> <td>W-15/136</td> <td>Sewn</td> <td>Ultrafresh</td> </tr> <tr> <td>Carrflex : CFX7EF</td> <td>120lbs/sq.in</td> <td>500 G/m2/24hr</td> <td>4</td> <td>W-4.5/110</td> <td>RF Welded & Sewn</td> <td>Fresche</td> </tr> </tbody> </table>				Moisture Resistance	Vapor Exchange	Stretch Way	Tear/Break Strength	Seams Type	Antimicrobial Treatment	ChemflexII:	150lbs/sq.in	194 G/m2/24hr	2	W-15/136	Sewn	Ultrafresh	Carrflex : CFX7EF	120lbs/sq.in	500 G/m2/24hr	4	W-4.5/110	RF Welded & Sewn	Fresche
	Moisture Resistance	Vapor Exchange	Stretch Way	Tear/Break Strength	Seams Type	Antimicrobial Treatment																		
ChemflexII:	150lbs/sq.in	194 G/m2/24hr	2	W-15/136	Sewn	Ultrafresh																		
Carrflex : CFX7EF	120lbs/sq.in	500 G/m2/24hr	4	W-4.5/110	RF Welded & Sewn	Fresche																		
Disinfection	<p>Antimicrobially treated to inhibit colonization. Wipe or rinse to clean and disinfect between uses. Compatible with bleach, hydrogen peroxide, chlorhexidine based disinfectant wipes or soap & water. Do not machine wash. Do not submerge. Discard if cover damaged or contamination is suspected. We recommend discarding after use in COVID-19 isolation.</p>																							
<p>For best results, place cushion <u>under</u> incontinence pads and <u>over</u> slide/boost sheets. Please see "How to Place" and "How to Clean" instructions for details</p>																								
Ordering	<p>www.JewellNursingSolutions.com, www.BedSoreRescue.com Or call (650) 294-8557 Bulk pricing available</p>																							
Guarantee & Warranty	<p>3 months warranty for materials & construction under normal wear and tear (Excludes sharp objects cuts and temperatures >200 degrees) 10-day Satisfaction guarantee or your money back when purchasing 2 or less.</p>																							
<p>Meets flammability requirements of electronic & appliance repair, home furnishing and thermal insulation technical bulletin 117-2013 Contains NO added flame-retardant chemicals</p>																								
<p>Jewell Nursing Solutions is a nurse owned and operated company dedicated to bringing effective pressure injury prevention products and services to the entire caregiving community. BedSore Rescue™ Cushion (formally The BackBone™) is proudly designed, patented, and trademarked by Jewell Nursing Solutions. A certified woman owned business celebrating 10 years of helping people stop Pressure Injuries.</p>																								
<p>Made in the USA by certified manufacturer reg #CA39930 (IL). Distributed by Jewell Nursing Solutions</p>																								